

Torex Pro Ankle

- Removable malleolar shaped gel inserts cradle soft tissue structures •
- Adjustable elastic strapping allows dynamic compression •
- Boot provides strength and support •
- Reusable and long lasting •
- Non-toxic/latex-free •



APPLICATIONS

- Ankle sprain
- Swollen ankle
- Ankle surgery
- Arthritis
- Edema
- General ankle injuries requiring ice or heat

PRODUCT FEATURES

- Used Hot, MC2 provides soothing long lasting relief. Can be used repeatedly, saving hundreds of dollars over expensive one-time-use hot packs.
- Used Cold, MC2 provides fast first aid for acute soft tissue injuries and other conditions that require cold applications to combat swelling, pain and bruising. Store in the freezer or refrigerator.
- Includes insulating compression bandage to protect the skin from extreme temperature.

The Torex® family of innovative cold therapy products includes our patented radial roll-on sleeves which offer 360° coverage of hard to treat limbs and joints. Visit us at www.torexhealth.com to find out more.

Torex Pro Ankle

TXMC2
ONE SIZE
FITS MOST

torex®